

DNPH1 Rabbit pAb

货号: B20956

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,LO2,HeLa,Mouse heart,Mouse kidney,Rat kidney,Rat liver
细胞定位	cytosol,extracellular exosome,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-174 of human DNPH1 (NP_006434.1).
序列	MAAAMVPGRSESWERGEPRPALYFCGSIRGGREDRTLYERIVSRLRRFGTVLTEHVAAAELGARGEEAAGGDRLIHEQD LEWLQQADVVAEVTQPSLGVGYELGRAVAFNKRILCLFRPQSGRVLSAMIRGAADGSRFQVWDYEEGEVEALLDRYFEA DPPGQVAASPDPTT

靶点信息

研究背景	This gene was identified on the basis of its stimulation by c-Myc protein. The latter is a transcription factor that participates in the regulation of cell proliferation, differentiation, and apoptosis. The exact function of this gene is not known but studies in rat suggest a role in cellular proliferation and c-Myc-mediated transformation. Two alternative transcripts encoding different proteins have been described. [provided by Ref Seq, Jul 2008]
------	---

基因ID	10591
基因名	DNPH1
Swiss	O43598
别名	C6orf108;RCL;dj330M21.3;DNPH1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn