

ZNF460 Rabbit pAb

货号: B20960

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	64kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,Rat liver,Rat pancreas
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-210 of human ZNF460 (NP_006626.3).
序列	SEMQEYFLRPGTDPOSEKLHGKMSLEHEGLATADGICSMMIQNQVSPEDALYGFDSPVTDSLIHEGENSYKFEEMFNE NCFLVQHEQILPRVKPYDCPECGKAFGKSKH

靶点信息

研究背景	Zinc finger proteins, such as ZNF272, interact with nucleic acids and have diverse functions. The zinc finger domain is a conserved amino acid sequence motif containing 2 specifically positioned cysteines and 2 histidines that are involved in coordinating zinc. Kruppel-related proteins form 1 family of zinc finger proteins. See ZFP93 (MIM 604749) for additional information on zinc finger proteins.[supplied by OMIM, May 2004]
基因ID	10794

基因名	ZNF460
Swiss	Q14592
别名	HZF8;ZNF272;ZNF460

产品验证



Western blot analysis of ZNF460 expressed in BxPC-3,Rat liver,Rat pancreas using ZNF460 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn