

WBP4 Rabbit pAb

货号: B20972

产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat thymus
细胞定位	nuclear speck,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-376 of human WBP4 (NP_009118.1).
序列	MADYWKSQPKKFCDYCKCWIADNRPSVEFHERGKNHKENVAKRISEIKQKSLDKAKEEEKASKEFAAMEAAALKAYQEDL KRLGLESEILEPSITPVTSTIPTSTSNQQKEKKEKKRKKDPSKGRWVEGITSEGYHYYDLISGASQWEKPEGFQGDLKK TAVKTVWVEGLSEDGFTYYNTETGESRWEKPDDFIPHTSDLPSSKVNENSLGTLDESKSSDSHSDSDGEQAEEGGVS TETEKPKIKFKEKNKNSDGGSDPETQKEKSIQKQNSLGSNEEKSKTLKKSNPYGEWQEIKQEVESHEEVDLELPSTENEYV STSEADGGGEPKVVFEKTVTSLGVMADGVAPVFKRRRTENGKSRNLRQRGDDQ

靶点信息

研究背景	This gene encodes WW domain-containing binding protein 4. The WW domain represents a small and compact globular structure that interacts with proline-rich ligands. This encoded protein is a general spliceosomal protein that may play a role in cross-intron bridging of U1 and U2 snRNPs in the spliceosomal complex A. [provided by RefSeq, Jul 2008]
------	--

基因ID	11193
基因名	WBP4
Swiss	O75554
别名	FBP21;WBP4

产品验证



Western blot analysis of WBP4 expressed in Rat thymus using WBP4 Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn