

MLXIP Rabbit pAb

货号: B20976

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	91kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	mitochondrial outer membrane,nucleus
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human MLXIP (NP_055753.3).
序列	ARAHASAAPPPPRAGPGREEPPRRQIIHSGHFMVSSPHREHPPKKGDFDTVKNQTCQTYSFGKTSSCHLSIDASLTKLF ECMTLAYSGKLVSPKWKNFK

靶点信息

研究背景	This gene encodes a protein that functions as part of a heterodimer to activate transcription. The encoded protein forms a heterodimer with Max-like protein X (MLX) and is involved in the regulation of genes in response to cellular glucose levels. [provided by RefSeq, Mar 2014]
基因ID	22877

基因名	MLXIP
Swiss	Q9HAP2
别名	MIR;MONDOA;bHLHe36;MLXIP

产品验证



Western blot analysis of MLXIP expressed in Rat lung using MLXIP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn