

ARHGEF26 Rabbit pAb

货号: B21001

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	48kDa/88kDa/97kDa
实测分子量	102kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Mouse brain,Rat lung
细胞定位	cytosol
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 592-871 of human ARH GEF26 (NP_056410.3).
序列	ICQKTPKDSPKYEVCKRALKEVSKLVRLCNEGARKMERTEMMYTINSQLEFKIKPFPLVSSSRWLVKRGELTAYVEDTVLFS RRTSKQQVYFFLFNDVLIITKKKSEESYNVNDYSLRDQLLVESCDNEELNSSPGKNSSTMLYSRQSSASHLFTLTVLSNHA NEKVEMLLGAETQSERARWITALGHSSGKPPADRTSLTQVEIVRSFTAKQPDELSLQVADVLLIYQRVSDGWYEGERLRD GERGWFPMECAKEITCQATIDKNVERMGRLLGLETNV

靶点信息

研究背景	This gene encodes a member of the Rho-guanine nucleotide exchange factor (Rho-GEF) family. These proteins regulate Rho GTPases by catalyzing the exchange of GDP for GTP. The encoded protein specifically activates RhoG and plays a role in the promotion of macropinocytosis. Underexpression of the encoded protein may be a predictive marker of chemoresistant disease. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Oct 2011]
基因ID	26084
基因名	ARHGEF26
Swiss	Q96DR7
别名	ARHGEF26;CSGEF;HMFN1864;SGEF

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn