

ABCA12 Rabbit pAb

货号: B21003

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	cytoplasm,cytosol,plasma membrane
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 2456-2595 of human ABCA12 (NP_775099.2).
序列	CTRLAIMVNGKFQCIQSLQHIKSRFGRGFTVKVHLKNNKVTMETLTKFMQQLHFPKTYLKDQHLMSLEYHVPVTAGGVANIFDLLETNKTALNITNFLVSQTTLEEVFINFAKDQKSYETADTSSQGSTISVDSQDDQMES

靶点信息

研究背景	The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intracellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, and White). This encoded protein is a member of the ABC1 subfamily, which is the only major ABC subfamily found exclusively in multicellular eukaryotes. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2008]
基因ID	26154

基因名	ABCA12
Swiss	Q86UK0
别名	ARCI4A;ARCI4B;ICR2B;LI2;ABCA12

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn