

FEZF2 Rabbit pAb

货号: B21031

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	49kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Rat brain
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human FEZF2 (NP_060478.3).
序列	MASSASLETMVPPACPRAGASPATSKTLAFSIERIMAKTSEPRAPFEPRPGALEADGSQGKKLLNLCSP LPCMIPLQPLGYEV PSKTLLSYSELWKSSLRAGGGGGGGGGGGGGGGGGGAPVCGASGLCKTNCGVCKAELGLAPSALPAGRVIKPQVINQAV GLPASGSLYYFNYLDSTAYPPSELLSGHLFPSGLLNAQAPAALAAHPKLFLENAKLAGLAADKFPH PAPYPHKERLPAPLEQ VLKENSALTAERGGVKGH SKLPGGSADGKPNFTCE

靶点信息

研究背景	Transcription repressor. Required for the specification of corticospinal motor neurons and other subcerebr al projection neurons. May play a role in layer and neuronal subtype-specific patterning of subcortical proj ections and axonal fasciculation. Controls the development of dendritic arborization and spines of large la yer V pyramidal neurons. May be involved in innate immunity (By similarity.
------	---

基因ID	55079
基因名	FEZF2
Swiss	Q8TBJ5
别名	FEZ;FEZL;FKSG36;TOF;ZFP312;ZNF312;FEZF2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn