

HEMGN Rabbit pAb

货号: B21037

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	55kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,22Rv1,BxPC-3,SP2/0,Mouse spleen
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HEMGN (NP_932095 .1).
序列	MDLGKDQSHLKHHQTPDPHQEENHSPEVIGTWSLRNRELLRKRKAEVHEKETSQWLFGEQKKRKQRTGKGNRRGRK RQNTTELKVEPQPQIEKEIVEKA

靶点信息

研究背景	Regulates the proliferation and differentiation of hematopoietic cells. Overexpression block the TPA-induc ed megakaryocytic differentiation in the K562 cell model. May also prevent cell apoptosis through the act ivation of the nuclear factor-kappa B (NF-kB).
基因ID	55363

基因名	HEMGN
Swiss	Q9BXL5
别名	CT155;EDAG;EDAG-1;NDR;HEMGN;hemogen

产品验证



Western blot analysis of HEMGN expressed in Mouse spleen,HepG2,22Rv1 using HEMGN Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn