

RAB20 Rabbit pAb

货号: B21042

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	26kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	Cytoplasmic side,Cytoplasmic vesicle,Golgi apparatus,Lipid-anchor,phagosome,phagosome membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 80-190 of human RAB20 (NP_060287.1).
序列	YDVNHRQSLVELEDRLFGLTDTASKDCLFAIVGNKVDLTEEGALAGQEKEECSPNMDAGDRVSPRAPKQVQLEDAVALY KKILKYKMLDEQDVPAAEQMCFETSAKTGYNV

靶点信息

研究背景	Plays a role in apical endocytosis/recycling. Plays a role in the maturation and acidification of phagosome s that engulf pathogens, such as S.aureus and M.tuberculosis. Plays a role in the fusion of phagosomes wi th lysosomes.
基因ID	55647

基因名	RAB20
Swiss	Q9NX57
别名	RAB20

产品验证



Western blot analysis of RAB20 expressed in Mouse heart using RAB20 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn