

ZNF426 Rabbit pAb

货号: B21073

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	63kDa
实测分子量	63kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	TE-1,A-431
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 150-270 of human ZNF426 (NP_077011.1).
序列	QCGEVFSEHSCLKTHVRTQSTGNTHDCNQYGKDFLTLCCKTSTGEKLSEFNQSEKIFSLTPNIVYQRTSTQEKSFECSHCGKSFINESYLQAHMRTHNGEKLYEWRNYGPGFIDSTSLSVL

靶点信息

研究背景	Kaposi's sarcoma-associated herpesvirus (KSHV) can be reactivated from latency by the viral protein RTA . The protein encoded by this gene is a zinc finger transcriptional repressor that interacts with RTA to modulate RTA-mediated reactivation of KSHV. While the encoded protein can repress KSHV reactivation, RTA can induce degradation of this protein through the ubiquitin-proteasome pathway to overcome the repression. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2015]
------	---

基因ID	79088
基因名	ZNF426
Swiss	Q9BUY5
别名	K-RBP;ZNF426

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn