

ST6GALNAC5 Rabbit pAb

货号: B21086

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse testis,Mouse brain,Rat testis
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 28-336 of human ST6GALNAC5 (NP_112227.1).
序列	SLGGQKERPPQQQQQQQQQQQASATGSSQPAAESSTQQRPGVPAGPRPLDGYLGVADHKPLKMHCRDCALVTSSG HLLHSRQGSQIDQTECVIRMNDAPTRGYGRDVGNRRTSLRVIAHSSIQRILNRHDLLNVSQGTVFIFWGPSSYMRRDGKG QVYNNLHLLSQVLPRLKAFMITRHKMLQFDELFKQETGKDRKISNTWLSTGWFTMTIALELCDRINVYGMVPPDFCRDPN HPSVPYHYEFPFGPDECTMYLSHERGRKGSHHRFITEKRVFKNWARTFNIHFFQPDWKPESLAINHPENKPVF

靶点信息

研究背景	The protein encoded by this gene is a Golgi type II transmembrane glycosyltransferase. The encoded protein catalyzes the transfer of sialic acid to cell surface proteins to modulate cell-cell interactions. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016]
------	---

基因ID	81849
基因名	ST6GALNAC5
Swiss	Q9BVH7
别名	SIAT7-E;SIAT7E;ST6GalNAcV

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn