

SPZ1 Rabbit pAb

货号: B21095

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	
实测分子量	50KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	DU145
细胞定位	cytoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-430 of human SPZ1 (NP_115956.3).
序列	MASSAKSAEMPTISKTVNPTDPHQEYLDPRITIALFEIGSHSPSSWGSPLFLKNSSHQVTEQQTAQKFNNLLKEIKDILKN MAGFE EKITEAKELFEETNITEDVSAHKENIRGLDKINEMLSTNLPVSLAPEKEDNEKKQEMILETNITEDVSAHKENIRGLDKI NEMLSTNLPVSLAPEKEDNEKKQQMIMENQNSENTAQVFARDLVNRLEEKKVLNETQQSQEKAKNRLNVQEETMKIRNN MEQLLQEAHWSKQHTELSKLIQSYQKSQKDISETLGNGVGFGTQPNNEVSAKHELEEQVKKLSHDTYSLQLMAALLEN ECQILQQRVEILKELHHQKQGTLEKPIQINYKQDKKNQKPSEAKKVEMYKQNKQAMKGTFWKKDRSCRSLDVCLNKKA CNTQFNIHVARKALRGKMRSASSLR

靶点信息

研究背景	This gene encodes a bHLH-zip transcription factor which functions in the mitogen-activate protein kinase (MAPK) signaling pathway. Because of its role in the upregulation of cell proliferation and tumorigenesis, t his gene may serve as a target for Ras-induced tumor treatments. [provided by RefSeq, Oct 2011]
------	---

基因ID	84654
基因名	SPZ1
Swiss	Q9BXG8
别名	NYD-TSP1;PPP1R148;SPZ1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn