

NEK9 Rabbit pAb

货号: B21104

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	107kDa
实测分子量	107kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Raji,Mouse heart,Rat thymus,Rat lung,Rat spleen
细胞定位	cytosol,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 730-979 of human NEK9 (NP_149107.4).
序列	SKTIRSNSSGLSIGTVFQSSSPGGGGGGGGGEEEDSQQESETPDPSSGGFRGTMEADRGMEGLISPTEAMGNSNGASSS CPGWLRKELENAEFIPMPDPSPLSAAFSESEKDTLPYEELQGLKVASEAPLEHKPQVEASSPRLNPAVTCAGKGTPLTTPPA CACSSLQVEVERLQGLVLKCLAEQQLQQENLQIFTQLQKLNKKLEGGQQVGMHSGTQTAKEEMEMDPKPDLDSDSW CLLGTDSCRPSL

靶点信息

研究背景	This gene encodes a member of the NimA (never in mitosis A) family of serine/threonine protein kinases. The encoded protein is activated in mitosis and, in turn, activates other family members during mitosis. T his protein also mediates cellular processes that are essential for interphase progression. [provided by Re fSeq, Jul 2016]
------	--

基因ID	91754
基因名	NEK9
Swiss	Q8TD19
别名	APUG;LCCS10;NC;NERCC;NERCC1;NEK9

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn