

CDHR1 Rabbit pAb

货号: B21114

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,Mouse brain,Mouse thymus,Rat brain,Rat thymus
细胞定位	photoreceptor outer segment membrane,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human CDHR1 (NP_001165442.1).
序列	QANFAPHFFDNGVGSTNGNMALFSLPEDTPVGSHVYTLNGTDPEGDPISYHISFDPSTRSVFSVDPTFGNITLVEELDRER EDEIEAIIISISDGLNLVAEKVVILVTDANDEAPRFIQEPYVALVPEDIPAGSIIFKVHAVDRDTGSGGSVTYFLQNLHSPFAVDR HSGVLRQLQAGATLDYERSRTHYITVVAKDGGGRLHGADVVFSAATTTVTNVEDVQDMAPVFGTPYYGYVYEDTLPGSE VLKVAMDGDRGKPNRILYSLVNGNDGAFEINETSG

靶点信息

研究背景	This gene belongs to the cadherin superfamily of calcium-dependent cell adhesion molecules. The encoded protein is a photoreceptor-specific cadherin that plays a role in outer segment disc morphogenesis. Mutations in this gene are associated with inherited retinal dystrophies. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2013]
------	--

基因ID	92211
基因名	CDHR1
Swiss	Q96JP9
别名	CORD15;PCDH21;PRCAD;RP65;CDHR1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn