

# STK35 Rabbit pAb

货号: B21122

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Mouse                                                                                                |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000                                                                            |
| 理论分子量 |                                                                                                      |
| 实测分子量 | 58kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Mouse kidney                                                                                         |
| 细胞定位  | cytoplasm,nuclear body,nucleolus,nucleoplasm,nucleus                                                 |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                                                                                                |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 190-380 of human STK35 (NP_543026.2).                                                                            |
| 序列   | RPEGGGGSARPRYSLLAEIGRGSYGVVYEAVAGRSGARVAVKKIRCDAPENVELALAEFWALTSLKRRHQNVVQFEECVLQRNGLAQRMESHGNKSSQLYLRLVETSLKGERILGYAEEPCYLWFMFCEGGDLNQYVLSRRPDPATNKSFMLQLTSAIAFLHKNHIVHRDLKPDNILITERSGTPILKVA |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                             |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is a kinase that is predominantly found in the nucleus. However, it can interact with PDLIM1/CLP-36 in the cytoplasm and localize to actin stress fibers. The encoded protein may be a regulator of actin stress fibers in nonmuscle cells. [provided by RefSeq, Jul 2008] |
| 基因ID | 140901                                                                                                                                                                                                                                                                                                      |

|       |                     |
|-------|---------------------|
| 基因名   | STK35               |
| Swiss | Q8TDR2              |
| 别名    | CLIK1;STK35L1;STK35 |

产品验证



Western blot analysis of STK35 expressed in Mouse kidney using STK35 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)