

EML3 Rabbit pAb

货号: B21136

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	95kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1
细胞定位	cytoplasm,microtubule cytoskeleton,nucleus,spindle,spindle microtubule
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 300-530 of human EML3 (NP_694997.2).
序列	HYRGHTDCVRCLAVHPDGVRVASGQTAGVDKDGKPLQPVVHIWDSETLLKLQEIGLGAFERGVGALAFSAADQGAFLCV VDDSNEHMLS VWDCSRGMKLAIEIKSTNDSVLAVGFNPRDSSCIVTSGKSHVHFWNWSGGVGVPNGNLTTRKQGVFGK YKKPKFIPCFVFLPDGDILTGDSEGNILTWGRSPSDSKTPGRGGAKETYGIVAQAHAHEGSIFALCLRRDGTVLS

靶点信息

研究背景	Regulates mitotic spindle assembly, microtubule (MT)-kinetochore attachment and chromosome separation n via recruitment of HAUS augmin-like complex and TUBG1 to the existing MTs and promoting MT-based MT nucleation. Required for proper alignment of chromosomes during metaphase.
基因ID	256364

基因名	EML3
Swiss	Q32P44
别名	ELP95;EML3

产品验证



Western blot analysis of EML3 expressed in THP-1 using EML3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn