

SNAI3 Rabbit pAb

货号: B21144

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	32kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse thymus,Rat thymus
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 7-186 of human SNAI3 (NP_840101.1).
序列	VKTHSSHRVPNYRRLETQREINGACSACGGLVVPLLPRDKEAPSVPGDLQPWDRSSAVACISLPLLPRIEEALGASGLDAL EVSEVDPRASAAIVPLKDSLNLHLNLPPLLVLPTRWSPTLGPDRHGAPEKLLGAERMPRAPGGFECFHCHKPYHTLAGLAR HRQLHCHLQVGRVFTCK

靶点信息

研究背景	SNAI3 is a member of the SNAIL gene family, named for the Drosophila snail gene, which plays roles in m esodermal formation during embryogenesis (Katoh and Katoh, 2003 [PubMed 12579345]).[supplied by O MIM, Apr 2009]
基因ID	333929

基因名	SNAI3
Swiss	Q3KNW1
别名	SMUC;SNAIL3;ZNF293;Zfp293;SNAI3

产品验证



Western blot analysis of SNAI3 expressed in A-549,Mouse thymus,Rat thymus using SNAI3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn