

LIN9 Rabbit pAb

货号: B21156

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	57kDa/61kDa/63kDa
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis
细胞定位	nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 350-470 of human LIN9 (NP_775106.2).
序列	TLGGFPVEFLIQVTRLSKILMIKKEHIKKLREMNTAEKLSYSMPISIEFQRRYATIVLELEQLNKDLNKLHKVQQYCYELAPDQGLQPADQPTDMRRRCEEEAEIVRHANSSTGQPC

靶点信息

研究背景	This gene encodes a tumor suppressor protein that inhibits DNA synthesis and oncogenic transformation through association with the retinoblastoma 1 protein. The encoded protein also interacts with a complex of other cell cycle regulators to repress cell cycle-dependent gene expression in non-dividing cells. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]
基因ID	286826

基因名	LIN9
Swiss	Q5TKA1
别名	BARA;BARPsv;Lin-9;TGS;TGS1;TGS2;LIN9

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn