

TBCCD1 Rabbit pAb

货号: B21194

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	
实测分子量	70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,U-87MG,Mouse brain,Mouse spleen
细胞定位	cytoplasm,spindle pole centrosome
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human TBCCD1 (NP_060608.1).
序列	LTPTRPLILSGNQTVTFAPFHTHYPMLDHMARTGLATVPNYWDNPMVVCRENSDTRVFQLLPCEFYVFIIPFEMEGDTTEIPGGLPSVYQKALGQREQK

靶点信息

研究背景	Plays a role in the regulation of centrosome and Golgi apparatus positioning, with consequences on cell shape and cell migration.
基因ID	55171

基因名	TBCCD1
Swiss	Q9NVR7
别名	TBCCD1

产品验证



Western blot analysis of TBCCD1 expressed in Mouse brain,HeLa,U-87MG using TBCCD1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn