

PHLPP2 Rabbit pAb

货号: B21312

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	147kDa
实测分子量	140kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse eye,Rat eye
细胞定位	cytoplasm,cytosol,intercellular bridge,mitotic spindle,nucleoplasm,photoreceptor outer segment membrane
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1250-1323 of human PHLPP2 (NP_055835.2).
序列	DSLNLIEVATEVPKRKTGYFAAPTQMEPEDQFVVPHDLEEEVKEQMKQHQDSRLEPEPHEEDRTEPPEEFDAL

靶点信息

研究背景	PHLPP2 (PH domain and leucine rich repeat protein phosphatase 2) is a protein-coding gene. GO annotations related to this gene include protein serine/threonine phosphatase activity. An important paralog of this gene is PHLPP1.
基因ID	23035

基因名	PHLPP2
Swiss	Q6ZVD8
别名	PHLPP2;PHLPPL;PPM3B

产品验证



Western blot analysis of PHLPP2 expressed in Mouse eye,Rat eye using PHLPP2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn