

BCAR1 Rabbit pAb

货号: B21347

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:100 - 1:500 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	93kDa
实测分子量	130KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse brain
细胞定位	Cell junction,Cytoplasm,focal adhesion
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 495-612 of human BCAR1 (NP_055382.2).
序列	EPQEPLVQDLQAAVAQSAVHELLEFARSAVGNAAHSTDRALHAKLSRQLQKMEDVHQTTLVAHGQALDAGRGGSGAT LEDLDRLVACSRVPEDAKQLASFLHGNASLLFRRTKATA

靶点信息

研究背景	The protein encoded by this gene is a member of the Crk-associated substrate (CAS) family of scaffold proteins, characterized by the presence of multiple protein-protein interaction domains and many serine and tyrosine phosphorylation sites. The encoded protein contains a Src-homology 3 (SH3) domain, a proline-rich domain, a substrate domain which contains 15 repeat of the YxxP consensus phosphorylation motif for Src family kinases, a serine-rich domain, and a bipartite Src-binding domain, which can bind both SH2 and SH3 domains. This adaptor protein functions in multiple cellular pathways, including in cell motility, apoptosis and cell cycle control. Dysregulation of this gene can have a wide range of effects, affecting different pathways, including cardiac development, vascular smooth muscle cells, liver and kidney function, endothelial migration, and cancer.
基因ID	9564
基因名	BCAR1
Swiss	P56945
别名	BCAR1;CAS;CAS1;CASS1;CRKAS;P130Cas

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn