

OSR2 Rabbit pAb

货号: B21353

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	35kDa
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-312 of human OSR2 (NP_001135934.1).
序列	MGSKALPAPIPLHPSLQLTNYSFLQAVNTFPATVDHLQGLYGLSAVQTMHMHWTLGYPNVHEITRSTITEMAAAQGLVD ARFPFPALPFTTHLFHPKQGAIAHVLPALHKDRPRDFANLAVAATQEDPPKMGDLSKLSPLGSPISGLSKLTPDRKPSR GRLPSKTKKEFICKFCGRHFTKSYNLLIHERHTHTDERPYTCDICHKAFRRQDHLRDHRYIHSKEKPFKCCQECGKGFCQSRTL AVHKTLMQESPHKCPTCGRTFNQRSNLKTHLLTHTDIKPYSCEQCGKVFRRNCDLRRHSLTHTPRQDF

靶点信息

研究背景	OSR2 is a mammalian homolog of the Drosophila odd-skipped family of transcription factors (Lan et al., 2004 [PubMed 15175245]).
基因ID	116039

基因名	OSR2
Swiss	Q8N2R0
别名	OSR2

产品验证



Western blot analysis of OSR2 expressed in 293T using OSR2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn