

# FSTL4 Rabbit pAb

货号: B21373

产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                      |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IF/ICC                                                                                            |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200                                             |
| 理论分子量 | 93kDa                                                                                                |
| 实测分子量 | 85kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | U-87MG,U-251MG,SH-SY5Y,NCI-H460,A-549,Mouse heart,Mouse brain,Rat heart,Rat brain                    |
| 细胞定位  | extracellular region                                                                                 |
| 纯化    | Affinity purification                                                                                |

抗原信息


|      |                                                                                                                                                                                                                                            |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 22-250 of human FSTL4 (NP_055897.1).                                                                                                                         |
| 序列   | GWMDPGTSRGPDVGVGESQAEEPRSFVTRREGLSSHNELLASCGKKFCSRGSRCVLSRKTGEPECQCLEACRPSYVPV<br>CGSDGRFYENHCKLHRAACLLGKRITVIHSDKCFLKGDCTMAGYARLKNVLLALQTRLQPLQEGDSRQDPASQKRLLVE<br>SLFRDLADGNGHLSSSELAQHVLKKQDLDEDLLGCSPGDLLRFDDYNSDSSLTLREFYMAFQVVQLSLA |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | Predicted to enable brain-derived neurotrophic factor binding activity and calcium ion binding activity. Predicted to be involved in cell differentiation. Predicted to act upstream of or within negative regulation of brain-derived neurotrophic factor receptor signaling pathway; negative regulation of collateral sprouting; and negative regulation of dendritic spine development. Predicted to be located in extracellular region and secretory granule. |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |        |
|-------|--------|
| 基因ID  | 23105  |
| 基因名   | FSTL4  |
| Swiss | Q6MZW2 |
| 别名    | FSTL4  |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)