

DIP2C Rabbit pAb

货号: B21462

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	170kDa
实测分子量	155kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Rat lung,Rat heart,Rat brain
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-260 of human DIP2 C (NP_055789.1).
序列	LAKHKERKMAVPMPKRRSLVVQTSMDAYTPPDTSSGSEDEGSVQGDSQGTPTSSQGSINMEHWISQAIHGSTTSTSS SSTQSGGSGAAHRLADVMAQTHIENHSAPPDVTITYTSEHSIQVERPQGSGSRTAPKYGNAELMETGDGVPVSSRVSAK IQQ

靶点信息

研究背景	This gene encodes a member of the disco-interacting protein homolog 2 family. The protein shares strong similarity with a Drosophila protein which interacts with the transcription factor disco and is expressed in the nervous system.
基因ID	22982

基因名	DIP2C
Swiss	Q9Y2E4
别名	KIAA0934;DIP2C

产品验证



Western blot analysis of DIP2C expressed in Rat lung,293T,Rat heart using DIP2C Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn