

RBMY1A1 Rabbit pAb

货号: B21612

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC IF/ICC
推荐浓度	IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	nuclear speck,nucleolus,nucleoplasm
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 90-300 of human RBMY 1A1 (NP_005049.1).
序列	SGGRRRPPASSRNRSPSGSLRSARGSRGGTRGWLPSHEGHLDDGGYTPDLKMSYSRGLIPVKRGPSSRSGGPPPKKSAPSAVARNSWMGSQGPMsQRRENYGVPPRRATISSWRNDRMSTRHDGYATNDGNHPSCQETRDYAPPSRGYAYRDNGHSNRDEHSSRGYRNRHSSRETRDYAPPSRGHAYRDYGHsSRDESYSRGYRNRSSR

靶点信息

研究背景	This gene encodes a protein containing an RNA-binding motif in the N-terminus and four SRGY (serine, arginine, glycine, tyrosine) boxes in the C-terminus. This protein is thought to function as a splicing regulator during spermatogenesis. Multiple closely related paralogs of this gene are found in a gene cluster in the AZFb azoospermia factor region of chromosome Y. Most of these related copies are thought to be pseudogenes, though several likely encode functional proteins. [provided by RefSeq, Mar 2016]
------	---

基因ID	5940
基因名	
Swiss	P0DJD3
别名	RBM;RBM1;RBM2;RBMγ;RBMγ1C;YRRM1;YRRM2;RBMγ1A1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn