

MKNK1 Rabbit pAb

货号: B21636

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	39kDa/47kDa/51kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,HepG2,Mouse pancreas,Rat skeletal muscle
细胞定位	Cytoplasm,Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 374-465 of human MKN K1 (NP_003675.2).
序列	VQGQAPEKGLPTPQVLQRNSSTMDLTLFAAEIALNRQLSQHEENELAEPEALADGLCSMKLSPPCKSRLARRRALAQAGRGEDRSPPTAL

靶点信息

研究背景	This gene encodes a Ser/Thr protein kinase that interacts with, and is activated by ERK1 and p38 mitogen-activated protein kinases, and thus may play a role in the response to environmental stress and cytokines. This kinase may also regulate transcription by phosphorylating eIF4E via interaction with the C-terminal region of eIF4G. Alternatively spliced transcript variants have been noted for this gene.
------	--

基因ID	8569
基因名	MKNK1
Swiss	Q9BUB5
别名	MNK1;MKNK1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn