

ANKRD46 Rabbit pAb

货号: B21664

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Mouse brain,Rat lung
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-232 of human ANKRD 46 (NP_001257308.1).
序列	MSYVFVNDSSQTNVPLLQACIDGDFNYSKRLLSESGFDPNIRDSRGRTGLHLAAARGNVDICQLLHKFGADLLATDYQGNT ALHLCGHVDTIQFLVSNGLKIDICNHQGATPLVLAKRRGVNKDVIRLLESLEEQEVKGFNRGTHSKLETMQTAESAMES HLLNPNLQQGEGVLSSFRTTWQEFVEDLGFWRVLLLIFFVIALLSLGIAYYRRTLRLGSFARQDRSRIQAI

靶点信息

研究背景	This gene encodes a protein containing multiple ankyrin repeats. Ankyrin domains function in protein-protein interactions in a variety of cellular processes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2012]
基因ID	157567

基因名	ANKRD46
Swiss	Q86W74
别名	ANK-S;GENX-115279;ANKRD46

产品验证



Western blot analysis of ANKRD46 expressed in Mouse lung,Mouse brain,Rat lung using ANKRD46 Rab bit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 3 0ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure tim e: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn