

TAZ1 Rabbit pAb

货号: B21670

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	191KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 13-276 of mouse TAZ1 (NP_001186665.1).
序列	EPGCKGTTTRGLRKRLPLTPEPGEGGGCSLEARGCEEESRQKQRMVTQASGREETEGDKLAKENGKITEAPSDDPQPGT DLVRKTSITSSESLQTVECSEFQNMFLQSLDKEELGEGIKRRMRIKKCKSLENPPLEITKNEATQNIKVEFQDELFKNTPKH SCNSLSPGVEKNCSFELHDYSFLHSEGCNNENNFEKPHDVCLHTEENSLKLKKENLRNLAEKDDTRTTTKLLKTEKSVIA SKLLEESHLYQKNNGLTSC

靶点信息

研究背景	Important for normal spermatogenesis and male fertility. Specifically required for progression to the post-meiotic stages of spermatocyte development. Seems to be necessary for normal expression levels of a number of testis-expressed gene transcripts, although its role in this process is unclear.
基因ID	375337

基因名	TOPAZ1
Swiss	Q8N9V7
别名	C3orf77;TAZ1

产品验证



Western blot analysis of TAZ1 expressed in Mouse brain using TAZ1 Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn