

MARK3 Rabbit pAb

货号: B22355

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	
实测分子量	80kD
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	human brain tissue, PC-3 cells
细胞定位	cytoplasm,cytosol,dendrite,extracellular exosome,plasma membrane
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MARK3 (NP_001122390.2).
序列	MSTRPLPTVNERTDENHTSHGDGRQEVTSTRSGARCRNSIASCADEQPHIGNYRLKTIKGKGNFAKVKLARHILTGREVAIKIIDKTQLNPTSLQKL

靶点信息

研究背景	The protein encoded by this gene is activated by phosphorylation and in turn is involved in the phosphorylation of tau proteins MAP2 and MAP4. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]
基因ID	4140

基因名	MARK3
Swiss	P27448
别名	CTAK1;KP78;PAR1A;Par-1a;MARK3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn