

# HOXB2 Rabbit pAb

货号: B22365

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	38kDa
实测分子量	35KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 50-140 of human HOXB 2 (NP_002136.1).
序列	EQTFPSLQPGASTLQRPRSQKRAEDGPALPPPPPPPLAAPPAPPEFPWMKEKKSAKKPSQSATSPSPAASAVPASGVGSP ADGLGLPEAGG

靶点信息

研究背景	This gene is a member of the Antp homeobox family and encodes a nuclear protein with a homeobox DN A-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded protein functions as a sequence-specific transcription factor that is involved in development. Increased expression of this gene is associated with pancreatic cancer.
基因ID	3212

基因名	HOXB2
Swiss	P14652
别名	K8; HOX2; HOX2H; Hox-2.8

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)