

HJURP Rabbit pAb

货号: B22375

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	84kDa
实测分子量	80KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Mouse brain,Mouse spleen,Rat testis,Rat thymus
细胞定位	mitochondrion,nucleolus,nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human HJURP (NP_060880.3).
序列	MLGTLRAMEGEDVEDDQLLQKLRRSRRRFQRRMQRLIEKYNQPFEDTPVVQMATLTETPQGLRIWGGRLIKERNEGEIQ DSSMKPADRTDGSVQAAAWGPELPSHRTVLGADSKSGEVDATSDQEEVAVWALAPAVPQSPLKNELRRKYLTQVDILL QGAEYFECAGNRAGRDRVTPLPSLSPAVPA

靶点信息

研究背景	Centromeric protein that plays a central role in the incorporation and maintenance of histone H3-like variant CENPA at centromeres. Acts as a specific chaperone for CENPA and is required for the incorporation of newly synthesized CENPA molecules into nucleosomes at replicated centromeres. Prevents CENPA-H4 tetramerization and prevents premature DNA binding by the CENPA-H4 tetramer. Directly binds Holliday junctions.
------	---

基因ID	55355
基因名	HJURP
Swiss	Q8NCD3
别名	FAKTS; URLC9; hFLEG1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn