

Btbd11 Rabbit pAb

货号: B22402

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	121kDa
实测分子量	121KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of mouse Btbd11 (NP_082985.2).
序列	MARRGKKPVVRTLEDLTLD SGYGG AADSVRSSNLSLCCSDSHPASPYGGSCWPPLADSMHSRHNSFDTVNTALVEDSEGLDCAGQHCSRLLPDLDEVPWTLQELELLLLRSRDPRAGPAV

靶点信息

研究背景	Predicted to enable protein heterodimerization activity. Acts upstream of or within SMAD protein signal transduction. Predicted to be located in membrane. Predicted to be integral component of membrane. Is expressed in central nervous system; dorsal root ganglion; genitourinary system; and thymus primordium. Orthologous to human BTBD11 (BTB domain containing 11).
基因ID	74007

基因名	
Swiss	Q6GQW0
别名	6330404E16Rik

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn