

FXYD3 Rabbit pAb

货号: B22633

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	IHC: 1:50 - 1:200
理论分子量	9kDa
实测分子量	9kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	extracellular exosome,plasma membrane
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-87 of human FXYD3 (NP_005962.1).
序列	MQKVTLGLLVFLAGFPVLDANDLEDKNSPFYYDWHSLQVGG LICAGVLCAMGIIIVMSAKCKCKFGQKSGHHPGETPPLIT PGSAQS

靶点信息

研究背景	This gene belongs to a small family of FXYD-domain containing regulators of Na ⁺ /K ⁺ ATPases which share a 35-amino acid signature sequence domain, beginning with the sequence PFXYD, and containing 7 invariant and 6 highly conserved amino acids. This gene encodes a cell membrane protein that may regulate the function of ion-pumps and ion-channels. This gene may also play a role in tumor progression. Alternative splicing results in multiple transcript variants encoding distinct isoforms.
基因ID	5349

基因名	FXYD3
Swiss	Q14802
别名	MAT8; PLML

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn