

# Phospho-NF2-S518 Rabbit pAb

货号: B22636

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	19-30kDa/59-72kDa
实测分子量	70KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cell projection,Cytoplasm,Cytoplasmic granule,Cytoplasmic side,Nucleus,Nucleus,Peripheral membrane protein,cytoskeleton,filopodium membrane,perinuclear region,ruffle membrane
纯化	Affinity purification

抗原信息


抗原信息	A synthetic phosphorylated peptide around S518 of human NF2 (NP_000259.1).
序列	RLSME

靶点信息

研究背景	<p>This gene encodes a protein that is similar to some members of the ERM (ezrin, radixin, moesin) family of proteins that are thought to link cytoskeletal components with proteins in the cell membrane. This gene product has been shown to interact with cell-surface proteins, proteins involved in cytoskeletal dynamics and proteins involved in regulating ion transport. This gene is expressed at high levels during embryonic development; in adults, significant expression is found in Schwann cells, meningeal cells, lens and nerve. Mutations in this gene are associated with neurofibromatosis type II which is characterized by nervous system and skin tumors and ocular abnormalities. Two predominant isoforms and a number of minor isoforms are produced by alternatively spliced transcripts.</p>
------	--

基因ID	4771
基因名	NF2
Swiss	P35240
别名	NF2;ACN;BANF;SCH;merlin

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)