

PKMYT1 Rabbit pAb

货号: B22671

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	54kDa
实测分子量	54KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse testis,Mouse thymus,Rat testis
细胞定位	cytosol,endoplasmic reticulum,Golgi apparatus,nucleolus,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 75-362 of human PKMYT 1 (NP_004194.3).
序列	HQLQPRRVSVFRGEASETLQSPGYDPSRPESFFQQSFQRLSRLGHGSYGEVFKVRSKEDGRLYAVKRSMSPFRGPKDRAR KLAEVGSHEKVGQHPCCVRLEQAWEEGGILYLQTELCGPSLQQHCEAWGASLPEAQVWGYLRDTLLALAHLSQGLVHL DVKPANIFLGPRGRCKLGDFLLVELGTAGAGEVQEGDPYMAPELLQGSYGTAAADVFSLGLTILEVACNMELPHGGEG WQQLRQGYLPPEFTAGLSSELRSVLVMMLEPDKLRATAEALLALPVLRQP

靶点信息

研究背景	This gene encodes a member of the serine/threonine protein kinase family. The encoded protein is a membrane-associated kinase that negatively regulates the G2/M transition of the cell cycle by phosphorylating and inactivating cyclin-dependent kinase 1. The activity of the encoded protein is regulated by polo-like kinase 1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	---

基因ID	9088
基因名	PKMYT1
Swiss	Q99640
别名	MYT1; PPP1R126

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn