

NEURL1B Rabbit pAb

货号: B22673

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	59kDa
实测分子量	59KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	actin cytoskeleton,cytosol,early endosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 270-440 of human NEURL1B (NP_001136123.1).
序列	ASSPALLEADLRFHATRGPDVLSADRKVACAPRPDGGRTLVSERPLRPGESLFVEVGRPGLAAPGALAFGITSCDPGVL RPNELPADPDALLDRKEYWVVARAGVPVPSGGDALSFTRLPGGDVLLGINGRPRGRLLCVDTTQALWAFFAVRGGVAGQL RLLGTLQSSPA

靶点信息

研究背景	E3 ubiquitin-protein ligase involved in regulation of the Notch pathway through influencing the stability and activity of several Notch ligands.
基因ID	54492

基因名	NEURL1B
Swiss	A8MQ27
别名	neur2; NEURL3; RNF67B; hNeur2

产品验证



Western blot analysis of NEURL1B expressed in Mouse brain using NEURL1B Rabbit pAb at 1:1000. Second ary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn