

AHNAK Rabbit pAb

货号: B22675

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC IF/ICC
推荐浓度	IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	actin cytoskeleton,cytoplasm,cytosol,extracellular exosome,focal adhesion,nucleus,plasma membrane,plasma membrane protein complex,T-tubule
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-149 of human AHNAK (NP_001611.1).
序列	MEKEETTRELLLPNWQGSGLTIAQRDDGVFVQEVTONSPAARTGVVKEGDQIVGATIFYFDNLQSGEVTQLLNTMGHHTVGLKLHRKGDRSPEPGQTWTREVFSSCSSEVVLSGDDEEYQRIYTTKIKPRLKSEDGVEGDLGETQSR

靶点信息

研究背景	The protein encoded by this gene is a large (700 kDa) structural scaffold protein consisting of a central domain with 128 aa repeats. The encoded protein may play a role in such diverse processes as blood-brain barrier formation, cell structure and migration, cardiac calcium channel regulation, and tumor metastasis. A much shorter variant encoding a 17 kDa isoform exists for this gene, and the shorter isoform initiates a feedback loop that regulates alternative splicing of this gene.
------	--

基因ID	79026
基因名	AHNAK
Swiss	Q09666
别名	

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn