

CMIP Rabbit pAb

货号: B22698

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	86kDa
实测分子量	70KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,Mouse testis,Mouse kidney,Rat kidney
细胞定位	cytosol,nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 640-773 of human CMIP (NP_938204.2).
序列	RLTLPSKSTDADLARLLSSGSFGNLENLSLAFTNVTSACAHEHLIKLP SLKQLNLWSTQFGDAGLRLLSEHLTMLQVLNLCET PVT DAGLLALSSMKSLCSLNMNSTKLSADTYEDLKAKLPNLKEVDVRYTEAW

靶点信息

研究背景	This gene encodes a c-Maf inducing protein that plays a role in T-cell signaling pathway. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
基因ID	80790

基因名	CMIP
Swiss	Q8IY22
别名	TCMIP

产品验证



Western blot analysis of CMIP expressed in Mouse testis,Rat kidney,K-562 using CMIP Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. B locking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn