

DLGAP3 Rabbit pAb

货号: B22709

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	106kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	cholinergic synapse,glutamatergic synapse,neuromuscular junction,plasma membrane,postsynaptic density,postsynaptic specialization
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DLGAP3 (NP_001073887.1).
序列	MRGYHGDRGSHPRPARFADQQHMDVGPAARAPYLLGSREAFSTEPRFCAPRAGLGHISPEGPLSLSEGPSVGPEGGPAGAGVGGSSTFPRMYPGQGPFD

靶点信息

研究背景	May play a role in the molecular organization of synapses and neuronal cell signaling. Could be an adapter protein linking ion channel to the subsynaptic cytoskeleton. May induce enrichment of PSD-95/SAP90 at the plasma membrane.
基因ID	58512

基因名	DLGAP3
Swiss	O95886
别名	DAP3; SAPAP3

产品验证



Western blot analysis of DLGAP3 expressed in Mouse brain using DLGAP3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn