

NUP88 Rabbit pAb

货号: B22785

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	83kDa
实测分子量	88KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCI-H460,Jurkat,U-251MG,Rat testis
细胞定位	cytosol,nuclear envelope,nuclear pore,nucleoplasm
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 502-741 of human NUP 88 (NP_002523.2).
序列	SPPLLCTREDVEVAESPLRVLAETPDSFEKHIRSILQRSVANPAFLKASEKDIAPPPPEECLQLLSRATQVFREQYILKQDLAKEE IQRRVKLLCDQKKKQLEDLSYCREERKSLREMAERLADKYEEAKEKQEDIMNRMKKLLHSFHSelpVLSDSERDMKKELQL IPDQLRHLGNAIKQVTMKKDYQQQKMEKVLSPKPTIILSAYQRKCIQSILKEEGEHIREMVKQINDIRNHVNF

靶点信息

研究背景	The nuclear pore complex is a massive structure that extends across the nuclear envelope, forming a gateway that regulates the flow of macromolecules between the nucleus and the cytoplasm. Nucleoporins, a family of 50 to 100 proteins, are the main components of the nuclear pore complex in eukaryotic cells. The protein encoded by this gene belongs to the nucleoporin family and is associated with the oncogenic nucleoporin CAN/Nup214 in a dynamic subcomplex. This protein is also overexpressed in a large number of malignant neoplasms and precancerous dysplasias. Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	4927
基因名	NUP88
Swiss	Q99567
别名	FADS4

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn