

Phospho-RIP3-S227 Rabbit pAb

货号: B22858

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	IHC: 1:50 - 1:200
理论分子量	25kDa/27kDa/56kDa
实测分子量	25kDa/27kDa/56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane,Cytoplasm,Mitochondrion,cytosol
纯化	Affinity purification

抗原信息

抗原信息	A synthetic phosphorylated peptide around S227 of human RIP3 (NP_006862.2).
序列	EPSLV

靶点信息

研究背景	The product of this gene is a member of the receptor-interacting protein (RIP) family of serine/threonine p rotein kinases, and contains a C-terminal domain unique from other RIP family members. The encoded pr otein is predominantly localized to the cytoplasm, and can undergo nucleocytoplasmic shuttling depende nt on novel nuclear localization and export signals. It is a component of the tumor necrosis factor (TNF) r eceptor-I signaling complex, and can induce apoptosis and weakly activate the NF-kappaB transcription fa ctor.
基因ID	11035

基因名	RIPK3
Swiss	Q9Y572
别名	RIPK3;RIP3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn