

CTIP2/BCL11B Rabbit pAb

货号: B22872

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	95kDa
实测分子量	120-130KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse thymus
细胞定位	neuron projection,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 190-310 of human CTIP 2/BCL11B (NP_612808.1).
序列	PVSGDGTQGEGQTEAPFGCQCQLSGKDEPSSYICTTCKQPFNSAWFLLQHAQNTHGFRIYLEPGPASSSLTPRLTIPPLGPEAVAQSPLMNFLGDSNPFNLLRMTGPILRDHPGFGEGR

靶点信息

研究背景	This gene encodes a C2H2-type zinc finger protein and is closely related to BCL11A, a gene whose translocation may be associated with B-cell malignancies. Although the specific function of this gene has not been determined, the encoded protein is known to be a transcriptional repressor, and is regulated by the NURD nucleosome remodeling and histone deacetylase complex. Four alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
------	--

基因ID	64919
基因名	BCL11B
Swiss	Q9C0K0
别名	ATL1; RIT1; CTIP2; IMD49; CTIP-2; IDDFSTA; ZNF856B; ATL1-beta; ATL1-alpha; ATL1-delta; ATL1-gamma; hRIT1-alpha

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn