

# SHANK3 Rabbit pAb

货号: B22875

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:1000 <b>IHC:</b> 1:50 - 1:200
理论分子量	184kDa
实测分子量	220KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	ciliary membrane,cytosol,dendritic spine,neuron projection,neuron spine,postsynaptic density
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human SHANK 3 (NP_277052.1).
序列	MDGPGASAVVVRVGIPDLQQTKCLRLDPAAPVWAAKQRVLCALNHSLQDALNYGLFQPPSRGRAGKFLDEERLLQEYPP NLDTPLPYLEFRYKRRVYAQNLIDDKQFAKLHTKANLKKFMDYVQLHSTDKVARLLDKGLDPNFHDPDSGE

靶点信息

研究背景	This gene is a member of the Shank gene family. Shank proteins are multidomain scaffold proteins of the postsynaptic density that connect neurotransmitter receptors, ion channels, and other membrane proteins to the actin cytoskeleton and G-protein-coupled signaling pathways. Shank proteins also play a role in synapse formation and dendritic spine maturation. Mutations in this gene are a cause of autism spectrum disorder (ASD), which is characterized by impairments in social interaction and communication, and restricted behavioral patterns and interests. Mutations in this gene also cause schizophrenia type 15, and are a major causative factor in the neurological symptoms of 22q13.3 deletion syndrome, which is also known as Phelan-McDermid syndrome. Additional isoforms have been described for this gene but they have not yet been experimentally verified.
基因ID	85358
基因名	SHANK3
Swiss	Q9BYB0
别名	PSAP2; SCZD15; PROSAP2; SPANK-2; DEL22q13.3

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)