

GTPBP10 Rabbit pAb

货号: B22929

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	42kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,HeLa,A-549
细胞定位	mitochondrion,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-387 of human GTPBP10 (NP_149098.2).
序列	MVHCSCVLFRKYGNFIDKLRLFTTRGGSGGMGYPRLGEGGGKGGDVWVVAQNRMTLKQLKDRYPRKRFVAGVGANSKISALKGSGKGDCEIPVPVGISVTDENGKIIGELNKENDRILVAQGGGLGGKLLTNFLPLKGQKRHHLDLKLADVGLVGFPNAGKSSLLSCVSHAKPAIADYAFTTLKPELGKIMYSDFKQISVADLPGLIEGAHMNKGGMGHKFLKHIERTRQLLFVVDISGFQLSSHHTQYRTAFETIILLTKELELYKEELQTKPALLAVNKMDLPDAQDKFHELMSQLQNPKDFLHLFEKNMIPERTVEFQHIIPISAVTGEIEELKNCIRKSLDEQANQENDALHKKQLNLWISDTMSSTEPPSKHAVTTSKMDII

靶点信息

研究背景	Small G proteins, such as GTPBP10, act as molecular switches that play crucial roles in the regulation of fundamental cellular processes such as protein synthesis, nuclear transport, membrane trafficking, and signal transduction
------	--

基因ID	85865
基因名	GTPBP10
Swiss	A4D1E9
别名	ObgH2; UG0751c10

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn