

HCoV-OC43 Nucleoprotein Rabbit pAb

货号: B22930

产品信息

反应	HCoV-OC43
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	49kDa
实测分子量	50KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human coronavirus (HCoV-OC43) Nucleocapsid Protein (His Tag)
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100-200 of coronavirus Nucleoprotein (YP_009555245.1).
序列	YWYRHNRRSFKTADGNQRQLLPWYFYLLGTGPHAKDQYGTIDIGVYWVASNQADVNTPADIVDRDPSSDEAIPTRFP PGTVLPQGYIEGSGRSAPNSRS

靶点信息

研究背景	Packages the positive strand viral genome RNA into a helical ribonucleocapsid (RNP) and plays a fundamental role during virion assembly through its interactions with the viral genome and membrane protein M. Plays an important role in enhancing the efficiency of subgenomic viral RNA transcription as well as viral replication.
基因ID	39105221

基因名	
Swiss	P33469
别名	

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn