

MERS-CoV Nucleoprotein Rabbit pAb

货号: B22931

产品信息

反应	MERS-CoV
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	44kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-413 of coronavirus Nu cleoprotein (YP_009047211.1).
序列	MASPAAPRAVSFADNNDITNTNLSRGRGRNPKPRAAPNNTVSWYTGLTQHGVPLTFPPGQGVPLNANSTPAQNAGYW RRQDRKINTGNGIKQLAPRWYFYTTGTGPEAALPFRVAVKDGIVVWHEDGATDAPSTFGTRNPNNDSAIVTQFAPGTKLPK NFHIEGTGGNSQSSSRASSLSRNSSRSSSQGSRSGNSTRGTSPGPSGIGAVGGDLLYDLLNRLQALESQKVKQSQPKVI TKKDAAAAKNMRHKRTSTKSFNMVQAFGLRGPGLQGNFGDLQLNKLGTEDPRWPQIAELAPTASAFMGMSQFKLTH QNNDDHGNPVYFLRYSGAIKLDPKNPYNNKWLELLEQNIDAYKTFPKKEKKQKAPKEESTDQMSEPPKEQRVQGSITQRT RTRPSVQPGPMIDVNTD

靶点信息

研究背景	Packages the positive strand viral genome RNA into a helical ribonucleocapsid (RNP and plays a fundame ntal role during virion assembly through its interactions with the viral genome and membrane protein M. Plays an important role in enhancing the efficiency of subgenomic viral RNA transcription as well as viral r eplication.
------	---

基因ID	14254601
基因名	
Swiss	K9N4V7
别名	

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn