

SARS-CoV-2 N Protein Rabbit pAb

货号: B22935

产品信息

反应	SARS-CoV-2
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	46kDa
实测分子量	55KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of coronavirus N Protein (NP_828858.1).
序列	MSDNGPQSNQRSAPRITFGGPTDSTDNNQNGGRNGARPKQRRPQGLPNNTASWFTALTQHGKEELRFPRGQGVPIINTNSGPDDQIGYYRRATRRVRGGDG

靶点信息

研究背景	Packages the positive strand viral genome RNA into a helical ribonucleocapsid (RNP and plays a fundame ntal role during virion assembly through its interactions with the viral genome and membrane protein M. Plays an important role in enhancing the efficiency of subgenomic viral RNA transcription as well as viral r eplication. May modulate transforming growth factor-beta signaling by binding host SMAD3.
基因ID	1489678

基因名	
Swiss	P59595
别名	

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn