

UGT2B17 Rabbit pAb

货号: B23064

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	61kDa
实测分子量	61kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 50-200 of human UGT2 B17 (NP_001068.1).
序列	GHEVIVLTSSASILVNASKSSAIKLEVYPTSLTKNDLEFFMKMFDRWTYSISKNTFWSYFSQLQELCWEYSDYNIKLCEDA VLNKKLMRKLQESKFDVLLADAVNPCGELLAELLNIPFLYSLRFSVGTYVEKNGGGFLFPPSYVPVMS

靶点信息

研究背景	This gene encodes a member of the uridine diphosphoglucuronosyltransferase protein family. The encoded enzyme catalyzes the transfer of glucuronic acid from uridine diphosphoglucuronic acid to a diverse array of substrates including steroid hormones and lipid-soluble drugs. This process, known as glucuronidation, is an intermediate step in the metabolism of steroids. Copy number variation in this gene is associated with susceptibility to osteoporosis.
基因ID	7367

基因名	UGT2B17
Swiss	O75795
别名	BMND12; UDPGT2B17

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn