

# FAM3A Rabbit pAb

货号: B23066

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                                |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB IHC                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200                                                |
| 理论分子量 | 21kDa/23kDa/25kDa                                                                                    |
| 实测分子量 | 21kDa/23kDa/25kDa                                                                                    |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  |                                                                                                      |
| 细胞定位  | Secreted                                                                                             |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                                                                                                              |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 34-230 of human FAM3 A (NP_068578.2).                                                                                          |
| 序列   | FPRIQQLFTSPESSVTAAPRARKYKGLPQPCPEEHLAFRVVSGAANVIGPKICLEDKMLMSSVKDENVGRGLNIALVNGVS<br>GELIEARAFDMWAGDVNDLLKFIRPLHEGTLVFVASYDDPATKMNEETRKLFSELGSRNAKELAFRDSWVFVGAKGVQN<br>KSPFEQHVKNKHSNKYEGWPEALEMEGCIPRRSTAS |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                         |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a cytokine-like protein. The expression of this gene may be regulated by peroxisome p roliferator-activated receptor gamma, and the encoded protein may be involved in the regulation of gluco se and lipid metabolism. Alternative splicing results in multiple transcript variants. |
| 基因ID | 60343                                                                                                                                                                                                                                                                                                   |

|       |                              |
|-------|------------------------------|
| 基因名   | FAM3A                        |
| Swiss | P98173                       |
| 别名    | FAM3A;2.19;DLD;DXS560S;XAP-7 |

产品验证



Western blot analysis of FAM3A expressed in human liver cancer using FAM3A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)