

Timm22 Rabbit pAb

货号: B23126

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:1000 IP: 1:100 - 1:500
理论分子量	20kDa
实测分子量	20KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	T47D,Mouse brain,Rat testis,Rat brain
细胞定位	mitochondrial inner membrane,mitochondrion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 19-121 of mouse Timm22 (NP_062792.2).
序列	AEAPLQYSLLLQYLVGDKRQPRLLPEGSLGGIPSPAKSEEQKMIERAMESCAFKAVLACVGGFVLGGAFGIFTAGIDTNVGFDPKDPYRTPTAKEVLKDMGQR

靶点信息

研究背景	Essential core component of the TIM22 complex, a complex that mediates the import and insertion of multi-pass transmembrane proteins into the mitochondrial inner membrane. In the TIM22 complex, it constitutes the voltage-activated and signal-gated channel. Forms a twin-pore translocase that uses the membrane potential as external driving force in 2 voltage-dependent steps (By similarity).
基因ID	56322

基因名	
Swiss	Q9CQ85
别名	Tim22

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn